

BLOOD BANK REFRIGERATOR BR-250



BLOOD BANK REFRIGERATOR available in three different sizes and designed to meet strict requirements for storage of whole blood. The blood selection of 4°C Refrigerators includes the Micro controlled temperature control having audio visual alarm in case of deviation.

*Extra-strength non CFC Refrigeration system.

*CFC free PUF(Poly Urethane Foam) insulation to ensure temperature stability and save power consumption.

Specification :

Model No. : BR 250

Temperature : 4°C

Accuracy : $\pm 0.5^\circ\text{C}$

Volume : 250 liter

No. of Baskets : Four

Blood Bags (Aprox) : 120

Dimension : 570 x 550 x 1420

Net Weight : 86 KG

Chart Recorder : Optional

Alarm : Audio / Visual

Light : Fluorescent

Operating Volt : 230V AC

Power Frequency : 50/60 Hz

Hi-Tech Instruments
Office 7 & 8 2nd Floor Block C-1,
Taj Medical Complex
M.A. Jinnah Road
Karachi Pakistan.
Phone: 0092 21 2787496
Mobile: 03002343299
E-Mail: hitech.aslam@gmail.com

