

- Blood Bank Refrigerator
- Blood Bag Shaker
- Tube Stripper
- Plasma Extractor
- Platelet Incubator
- Platelet Agitator
- Platelet Rotator
- Blood Donor Chair
- Plasma Thawing Bath
- Economy Incubator
- Economy Oven

BLOOD BANK REFRIGERATOR BR500

BLOOD BANK REFRIGERATOR available in six different sizes and designed to meet strict requirements for storage of whole blood. The blood selection of 4°C Refrigerators includes the Micro controlled temperature control having audio visual alarm in case of deviation.

*Extra-strength non CFC Refrigeration system.

*CFC free PUF(Poly Urethane Foam) insulation to ensure temperature stability and save power consumption.

Specification :

Model No. : BR 500
 Temperature : 4°C
 Accuracy : $\pm 0.5^\circ\text{C}$
 Volume : 500 liter
 No. of Baskets : Six
 Blood Bags (Aprox) : 250
 Dimension : 640 x 653 x 2050
 Net Weight : 96 KG
 Chart Recorder : Optional
 Alarm : Audio / Visual
 Light : Fluorescent
 Operating Volt : 230V AC
 Power Frequency : 50/60 Hz



Hi-Tech Instruments
 Office 7 & 8nd Floor Block C-1,
 Taj Medical Complex
 M.A. Jinnah Road
 Karachi Pakistan.
 Phone: 0092 21 2787496
 Mobile: 03002343299
 E-Mail: hitech.aslam@gmail.com