

## LABORATORY REFRIGERATOR LR-345



Laboratory refrigerator designed to meet strict requirement of the storage of vaccines, reagents and other critical samples. Protecting the integrity of your work is imperative which is why many research and clinical laboratories rely on us for our field of expertise in high quality refrigeration solutions. It provides superior temperature uniformity, reliability and an eco-friendly construction.

\*Micro controlled temperature control having audio visual alarm in case of deviation.

\*Extra-strength non CFC Refrigeration system.

\*CFC free PUF(Poly Urethane Foam) insulation to ensure temperature stability and save power consumption.

### Specification :

Model No. : LR 345

Temperature : 4°C

Accuracy :  $\pm 0.5^{\circ}\text{C}$

Volume : 345 liter

No. of Shelves : Five

Dimension : 756 x 685 x 2035

Net Weight : 86 KG

Alarm : Audio / Visual

Light : Fluorescent

Operating Volt : 230V AC

Power Frequency : 50/60 Hz

Hi-Tech Instruments  
Office 7 & 8 2<sup>nd</sup> Floor Block C-1,  
Taj Medical Complex  
M.A. Jinnah Road  
Karachi Pakistan.  
Phone: 0092 21 2787496  
Mobile: 03002343299  
E-Mail: hitech.aslam@gmail.com

