

PLATELET INCUBATOR

PLATELET INCUBATOR an ideal storing environment for the platelets. This unit is designed to maintain a stable and uniform temperature of 22.0°C.. Simple calibration is possible via a control panel this coupled with a visual and audible alarm system, ensures that your samples are both safe and secure.

Model No. PI 900
Operating Volt 230V AC
Power Frequency 50/60 Hz
Power Consumption 15 watts

PID Microcontroller Based Temperature Controller
Stainless Steel Interior
Powder Coated Exterior
Audio / Visual Alarm
Fluorescent Light
CFC free cooling System



Hi-Tech Instruments
Office 7 & 8th Floor Block C-1,
Taj Medical Complex
M.A. Jinnah Road
Karachi Pakistan.
Phone: 0092 21 2787496
Mobile: 03002343299
E-Mail: hitech.aslam@gmail.com

