

## PLATELET INCUBATOR

PLATELET INCUBATOR an ideal storing environment for the platelets. This unit is designed to maintain a stable and uniform temperature of 22.0°C. Simple calibration is possible via a control panel this coupled with a visual and audible alarm system, ensures that your samples are both safe and secure.

Model No. PT 950  
Operating Volt 230V AC  
Power Frequency 50/60 Hz  
Power Consumption 15 watts

PID Microcontroller Based Temperature Controller
Powder Coated Interior
Powder Coated Exterior
Audio / Visual Alarm
Fluorescent Light
CFC free cooling System



Hi-Tech Instruments  
Office 7 & 8<sup>th</sup> Floor Block C-1,  
Taj Medical Complex  
M.A. Jinnah Road  
Karachi Pakistan.  
Phone: 0092 21 2787496  
Mobile: 03002343299  
E-Mail: hitech.aslam@gmail.com

